



**"Better Hearing,  
Better Living."**

### Contact Us

 JSS Institute of Speech and Hearing  
MG Road, Mysuru – 570004

 [jssish@jssonline.org](mailto:jssish@jssonline.org)  
[jssishmys@rediffmail.com](mailto:jssishmys@rediffmail.com)

 0821-2548229, 0821-2335014



JSS Mahavidyapeetha

**JSS INSTITUTE OF SPEECH AND HEARING**

Mahatma Gandhi Road, Mysuru 570004

**HOW TO USE COCHLEAR IMPLANT**



# HOW TO USE COCHLEAR IMPLANT

- Press the on/off button down for 2 sec.
- Place the processor on your ear, letting the coil dangle. Move the coil sideways and onto your implant.
- To change the program, press and hold the program button down for 2 sec. The number of times the light blinks indicates the program number.
- Press the top volume button to turn up the volume or sensitivity and the lower button to turn it down.
- To change Coil cable Hold the grips and pull the coil cable straight out without twisting. Then, hold the plug and pull the coil cable out of the coil. To put in the new coil cable, push the plug into the coil until you hear a click. Then, push the coil cable into the processing until it clicks.
- Replace ear hook if it gets loose. To replace the ear hook, pull it up to remove and push the new one in.
- To Change Microphone guard every 3 months. Use your fingernail to remove the old cover. Clip the new one on with the narrow side facing the ear-hook.
- To remove the battery, twist the processor and battery at an angle. Take out the battery and replace it with a charged one.
- To change the magnet, twist it left to take it out. Put the new magnet in the hole from the top. The side with the grips should face up. Twist the magnet right until it clicks into place.

## Charging :

- To recharge the battery, twist the processor and battery to remove the module. Place the battery on the charger at an angle. Twist it right to connect it to the charger. Plug the power adaptor into the socket and turn it on. When the charger light turns green, the battery is fully charged."



## Accessories:



**Safety cord:** For safety, attach the safety line to the processor. Clip the other end to the child's clothes.



**Huggies:** Huggies will help processor stay on the child's ear.



**Baby Lite Wear:** Helps make the processor lighter for little kids to wear.



**Headbands:** Headbands can be used for retention of the processor.